



Heavy-Duty Rumen Magnet



Weiyou® Heavy-duty rumen magnet is made of high quality Ferrite Magnet. We have been specialized in veterinary products for many years. Weiyou® Heavy-duty rumen magnet have a low-price advantage and cover most of the European and American markets. Weiyou® is committed to supplying leading-edge livestock identification and animal husbandry products.

Product Description

Weiyou® Heavy-duty rumen magnet is made of high quality Ferrite Magnet. Rumen Magnet with plastic cage is mainly used for cattle and cow, and its durable structure boasts plated, solid steel end caps for corrosion prevention and a solid steel connector rod through the core to reduce breakage in shipping and handling. Heavy-duty rumen magnet Stronger Magnetization before assembly creates a consistent pulling force from each section along the length of the magnet and enhances its overall strength.

Specification

size	D 0.69" X 3.1"-S, D 0.73" X 3.1"-M
Material	Ferrite Magnet.
Usage	cattle
Shelf life	3-5 years

Features

Durable Its durable structure boasts plated, solid steel end caps for corrosion prevention and a solid steel connector rod through the core to reduce breakage in shipping and handling. No plastic parts involved.

Stronger Magnetization before assembly creates a consistent pulling force from each section along the length of the magnet and enhances its overall strength.

Superior Characteristically, These qualities make it easier for the magnet to travel to the stomach and stay there.

size: dia 3/4" X 3-1/8" long

packaging: 12pcs per neuter packaging.





















