



Round Hub Veterinary Needles

Weiyou® round hub veterinary needles is made of SUS304 stainless steel which is reusable with many advantages. Weiyou® manufactures and supplies round hub veterinary needles of various sizes, and we has a low-price advantage. The round hub veterinary needles is used for injections in various animals, and it's very very sharp needle with tri-beveled.

Product Description

Weiyou® is a professional manufacturer in China, and we supply round hub veterinary needles which is very high quality. The round hub veterinary needles is made of SUS304 stainless steel pipe and luer lock brass round hub. With Ultra-sharp, tri-bevelled and Sterileneedle, which is best suitable for veterinary use. Leak proof joint between Hub and cannula prevents cannula come out of Hub during injection. Weiyou® manufactures and supplies round hub veterinary needles of various sizes, and we has a low-price advantage. Various sizes are available and OEM is acceptable.

Specification

Product name	Veterinary needles
Size	14G-24G
Feature	Mechanical
Material	Stainless steel
Usage	Veterinary syringe

Features

- SUS304 stainless steel pipe
- Luer lock brass round hub
- Pipe&Hub rivet connection in five direction
- Very sharp needle with tri-beveled
- 1 dozen per plastic box, 1 gross per white box
- 50 grosses with one export carton
- small or big knurl for your choice







www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone: +86-13567833512
E-mail: dario@nbweiyou.com















